




New Arrival!



The GelVet Max provides increased safety for returning current flow to the generator, and acts as a fully waterproof and reusable pressure reduction pad. It removes the need for disposable sticky plates during electro-surgical procedures, removing any risk of patient harm due to pad-site burns.

The Pad is filled with a viscoelastic polymer, which effectively reduces patient pressure sores associated with direct pressure, shearing and friction. The GelVet Max is ideal for a wide variety of surgical procedures making it a key item for both safety and simplicity during surgery.

No Need to Shave


Using the GelVet Max eliminates the need to shave the patient or check for metal work, avoid scar tissue or bony prominences. Just plug in and go.

No Pad-Site Burns

Unintended patient burns are a real concern with standard sticky plates. The solution is the GelVet Max, the safest return electrode available.

Reusable

Stop buying consumable single use REM plates and save money by using the reusable GelVet Max.



SKU	Description	Price exc GST
AHD-003	GelVet Max Reusable Return Electrode Pad - 85 x 40cm	\$4,100.00
AHD-CC-03	Microphone Head and Single Pole Connector	\$62.00
AHD-CC-05	Flat Head and Single Pole Connector	\$62.00
AHD-CC-06	Flat Head and Dual Pole Connector	\$62.00

Compatibility

The GelVet Max is compatible with most electrosurgical units. The GelVet Max does not include a connector, select your preferred connector type depending on your electrosurgery unit at the time of purchase.



For the full compatibility chart scan here



AHD-003
Microphone Head,
Single Pole



AHD-005
Flat Head,
Single Pole



AHD-006
Flat Head,
Dual Pole

ESU Brand	Recommendation
Zatha	Flat Head, Dual Pole
Valleylab	Flat Head, Dual Pole
ERBE	Flat Head, Single Pole
ConMed	Flat Head, Single Pole
Eschmann	Microphone Head, Single Pole